

# FUJINON Digital Video Processor

**New**

## EPX-2500



**FUJINON**  
FUJIFILM

# FUJINON Digital Video Processor

## EPX-2500

### Digital Processor and 150W Xenon Light Source Combination Unit


Fujinon introduces a higher standard for electronic video endoscopy with the newly designed EPX-2500.

Setting a new criteria for diagnostic and therapeutic endoscopy, the EPX-2500 provides physicians with Fujinon's world-class imaging and the exceptional ergonomics of the newly designed 500 series endoscopes.

## Specifications

Digital Output	DVI (Digital Visual Interface) : 1024x768p
Analog Output	RGB(2) : SDTV (NTSC / PAL) Y/C(2) : SDTV (NTSC / PAL) Composite : SDTV (NTSC / PAL)
Color Adjustment	Black, Red, Green, Blue, R-Hue, Chroma; 9 settings
Detail	Hi, LO; 9 settings
Contrast (Gamma)	9 settings
BLV	Hi, Mid, Lo, Off
Picture in Picture	On, Off; Size: 1/4, 1/3
Auto Gain Control	Off, +3db, +6db
Iris	Average/Peak
Zoom	Electric zoom : x1.0 - x2.0; 0.05steps
Lamp rated value	Main lamp : 11.7V 150W Xenon lamp Emergency lamp : 12V 75W Halogen lamp
Brightness control	9 settings
Lamp cooling method	Forced air cooling
Air supply pump	Hi, Low, Off
Power	120V 60Hz 2.7A / 230V 50Hz 1.4A
Dimensions (W x H x D)	375 x 190 x 495mm
Weight	17.0kg

## Function



FUJIFILM Group

### FUJINON Corporation

1-324 UETAKE, KITA-KU, SAITAMA CITY, SAITAMA 331-9624, JAPAN  
TEL : 81-48-668-2153 FAX : 81-48-668-1570  
E-mail: sales@msv.fujinon.co.jp URL <http://www.fujinon.co.jp/>

### FUJINON INC.

10 HIGH POINT DRIVE, WAYNE, NJ 07470, U.S.A.  
TEL : 1-800-490-0661 FAX : 1-973-633-8818  
E-mail : med.info@fujinon.com URL <http://www.fujinonendoscopy.com/>

### FUJINON (EUROPE) GmbH

HALSKESTRASSE 4, 47877 WILLICH, GERMANY  
TEL : 49-2154-924-0 FAX : 49-2154-924-290  
E-mail : Fujinon@fujinon.de URL <http://www.fujinon.de/>

### FUJINON SINGAPORE PTE. LTD.

BLK 211 HENDERSON ROAD, #10-04 HENDERSON INDUSTRIAL PARK,  
SINGAPORE, 159552 TEL : 65-6276-4988 FAX : 65-6276-6911  
E-mail : fujinon@fujinon.com.sg URL <http://www.fujinon.com.sg/>

### FUJINON AUSTRALIA PTY. LTD.

UNIT-18, 52 HOLKER STREET, SILVERWATER N.S.W. 2128, AUSTRALIA  
TEL : 61-2-9748-2744 FAX : 61-2-9748-2428 E-mail : sales@fujinon.com.au

### SHANGHAI FUJINON MEDICAL EQUIPMENT CO.,LTD.

FLAT B, 6/FL, SHI MEI BLD, 79 RI JING ROAD,  
WAI GAO QIAO FREE TRADE ZONE, SHANGHAI, ZIP CODE: 200131 CHINA  
TEL : 86-21-5424-9262(Office) TEL:86-21-5866-2151(Service)  
FAX : 86-21-5866-2150 E-mail : fjing2001@hotmail.com